

**CentriVet Blood Glucose and
Ketone
Monitoring System**

**Bovine Blood – CentriVet Blood
Ketone Strip**

Accuracy Study Report

Nov 2015

1.0 Objective

The objective of the accuracy study is to demonstrate that the accuracy of the CentriVet Blood Glucose and Ketone Monitoring System – CentriVet Blood Ketone Strip with bovine blood samples, when compared to a laboratory reference measurement method, meets the accuracy acceptance criteria.

2.0 Methods

Bovine blood samples from lactating cows were obtained. For each sample, concentration of D-3-hydroxybutyrate (β -ketone) was measured with a reference method and CentriVet Blood Glucose and Ketone Monitoring System with CentriVet Blood Ketone strip for comparison. Additionally, bovine blood samples were further spiked with D-3-hydroxybutyrate (β -ketone) to obtain bovine blood samples with high concentration of D-3-hydroxybutyrate (β -ketone).

Study Site:

The clinical study was conducted at AZURE Institute, San Diego, California, USA.

Sample Type:

Venous bovine blood sample was drawn from the tail vein of the subject with a syringe and needle and injected into a blood collection tube with heparin anticoagulant.

Reference Method:

Procedures were followed according to instructions given in Randox D-3-Hydroxybutyrate (Ranbut) assay kit package insert to prepare samples. Measurements were taken using Genesys UV-Vis Spectrophotometer.

Number of Strip Lots and Meters:

3 strip lots and 6 meters were tested for the study. The meters reported plasma equivalent blood β -ketone concentration values.

Code Chips:

Bovine code chips are specific for bovine blood sample with CentriVet Blood Ketone Strips.

3.0 Acceptance Criteria


95 % of the measured values shall fall within either ± 0.3 mmol/L of the average measured values of the reference measurement at D-3-hydroxybutyrate (β -ketone) concentrations < 1.5 mmol/L or within $\pm 20\%$ of the average measured values of the

reference measurement at D-3-hydroxybutyrate (β -ketone) concentrations $\geq 1.5\text{mmol/L}$.

4.0 Results

Regression Analysis

CentriVet Blood Glucose and Ketone Meter Reading – CentriVet Blood Ketone Strip vs. Laboratory Reference Method, Bovine Blood Sample: All 3 Strip Lots Combined



N = 336

Slope = 1.0208

Intercept = 0.0005

$R^2 = 0.9855$

Data Table

System Accuracy Results of CentriVet Blood Glucose and Ketone Monitoring System – Blood Ketone Strip, Bovine Blood Sample, for All 3 Lots Combined:

| CentriVet GK Monitoring System - Ketone Strip, Bovine Blood Sample | | | |
|--|-------------------------|---------------------------------------|---------------------|
| System Accuracy Results for D-3-Hydroxybutyrate Concentration ≥ 1.5 mmol/L | | | |
| Within $\pm 5\%$ | Within $\pm 10\%$ | Within $\pm 15\%$ | Within $\pm 20\%$ |
| 66 / 162 (40.7%) | 116 / 162 (71.6%) | 150 / 162 (92.6%) | 160 / 162 (98.8%) |
| System Accuracy Results for D-3-Hydroxybutyrate Concentration < 1.5 mmol/L | | | |
| Within ± 0.1 mmol/L | Within ± 0.2 mmol/L | Within ± 0.3 mmol/L | |
| 135 / 174 (77.6%) | 161 / 174 (92.5%) | 171 / 174 (98.3%) | |
| System Accuracy Results for D-3-Hydroxybutyrate Concentration ≥ 1.5 mmol/L and < 1.5 mmol/L | | | |
| Within $\pm 15\%$ or ± 0.3 mmol/L | | Within $\pm 20\%$ or ± 0.3 mmol/L | |
| 321 / 336 (95.5%) | | 331 / 336 (98.5%) | |

5.0 Conclusion

The results showed that more than 95% of data points for all 3 lots of CentriVet Blood Ketone Strip, with CentriVet Blood Glucose and Ketone Monitoring System, for bovine blood sample were within $\pm 20\%$ versus laboratory reference values when D-3-hydroxybutyrate (β -ketone) concentration is ≥ 1.5 mmol/L, or within ± 0.3 mmol/L versus laboratory reference values when D-3-hydroxybutyrate concentration is < 1.5 mmol/L. The system accuracy results meet the acceptance criteria.